

# The Mars Contrast Sensitivity Test

## USER MANUAL

### Two versions

**Mars Letter Contrast Sensitivity Test**

REF MLCST



**Mars Numeral Contrast Sensitivity Test**

REF MNCST



**mars perceptrix**

[www.marsperceptrix.com](http://www.marsperceptrix.com)  
Rev 05-06-2024



The Mars Perceptrix Corporation  
40 Valley View Road  
Chapel Hill, New York, 10514 USA  
Wellkang Ltd, Enterprise Hub,  
NW Business Complex, 1 Blanghorne Rd,  
Berry BT48 8SE, Northern Ireland






## What is the Mars Contrast Sensitivity Test and Why Use It?

**The Mars Contrast Sensitivity Test**, available in letter and optotype\* versions, consists of a set of calibrated, printed charts for assessment of visual contrast perception, and a simple, accurate scoring system. The test provides clinicians, researchers and occupational testers information about the minimum visual *contrast* required to support object perception. Most clinicians and researchers now believe that tests of contrast sensitivity are superior to ordinary visual acuity tests in predicting performance on many visual tasks of everyday life including reading, driving, and recognizing faces. With the Mars test, items gradually decrease in contrast instead of size, and what is evaluated is the minimum contrast a person needs for valid object identification.

Contrast sensitivity testing, a relatively recent advance in clinical and occupational vision assessment, gauges a person's ability to perceive low visual contrasts. Low contrasts are subtle lightness variations in the imagery we see, and they abound in the world around us (see a few examples in the pictures on the right). Contrast sensitivity loss occurs with most age-related eye disorders, including cataracts, glaucoma, and macular disease, and to a lesser degree, with normal aging. In some cases, contrast sensitivity loss occurs without reduction in visual acuity.

Difficulty seeing low contrasts makes performance of many everyday tasks more difficult. Many of these tasks are important for personal safety, employment, and mobility, as well as quality of life activities such as recreational activities, and appreciation of art and natural beauty. For many patients with vision impairment, reduced ability to discern low contrasts is often the perceptual change they notice most. In such cases it is important for clinicians to assess the severity of the loss in relation to patient complaints, to better understand and communicate with them about their visual function experience. Measuring contrast sensitivity with the Mars test provides an objective way to assess worsening or improving contrast vision over time. It is also useful for explaining to patients and their families how their contrast vision has declined, as visual disorders and low vision have progressed—or how they will improve after an intervention like cataract extraction. The Mars test therefore is an essential component in the toolbox of anyone who wants a more complete picture of visual functioning.

Vision professionals of all types (e.g. ophthalmologists, optometrists, orthoptists, occupational and low vision therapists), researchers, and occupational testers use the Mars test to (1) ascertain sufficient contrast sensitivity for specific task performance, such as driver's licensure, (2) establish baseline contrast sensitivity prior to an intervention, such as cataract extraction, (3) assess functional improvement in clinical trials, (4) identify functional losses in low contrast perception, (5) monitor progression of eye conditions, and (6) illustrate the concepts of image contrast and contrast sensitivity to patients and others.



---


\* The term *optotype* refers to a symbol that is used as a stimulus in testing vision function. Letters and numerals are the most common, but not the only optotypes.



illumination. When examinees hold the chart they generally will spontaneously adjust the chart position during testing and will aim the illumination onto the region containing the items they are trying to identify.

**Viewing distance and correction:** Testing should be performed with the eyes undilated to ensure a valid refractive state. The examinee's viewing distance should be at or about the customary near refraction distance of 40 cm (15.75 inches), but the test is tolerant of small refractive errors since the letters are large (20/600 equivalent at 40 cm). Examinees should wear their customary reading glasses. Alternatively, they can wear their distance correction supplemented with an additional +2.50 D, using a trial frame or Halberg clip. When testing monocularly, the untested eye should be covered with an occluder, patch, or hand. For patients with very low visual acuity who cannot easily read the highest contrast letters, test distance may be shortened to 20 cm (8 inches), increasing the add, if necessary, to +5.00 D; in this case care must be taken not to allow the patient's head to occlude the light source illuminating the chart.

### Testing Instructions



**Important:** If the examinee responds with an item other than one of those allowed, do not score the response as incorrect. Instead, remind the examinee of the allowed items, and ask for another response. This is in order to support the assumption that the probability of a guess is 1/10.

**Encourage the examinee to guess even when they report that the letters appear too faint.** This is to ameliorate the tendency of some subjects to “give up” before they have reported all the letters or numerals they are actually able to identify, or to be reluctant about giving incorrect answers when they are unsure of what they see. It helps ensure that the score is based on what the examinee *can* see and not on what the examinee *believes* he or she can see. Do not terminate the test because the examinee has given up and has stopped responding. Require the examinee to continue responding, even with guesses. Many times, such guesses reveal sensitivity to items beyond what the examinee believes they can identify. **Terminate testing only when the examinee makes two consecutive errors or reaches the end of the chart.**

## How to score the test

Scoring is simple. The test stops when the examinee has made two consecutive errors in identification. The log contrast sensitivity (**log CS**) score is simply the value accompanying the **final correct item** on the score sheet, minus a **scoring correction**. (If the examinee reaches the end of the chart without making two consecutive errors, then the final correct item is simply the last letter correctly identified.) The scoring correction is simply a penalty of 0.04 for each item missed prior to the final correct letter.

In the example below, there were no consecutive errors until the examinee failed to correctly identify the **O** at the end of line 6 and the **H** at the beginning of line 7. Thus, the final correct item is the **H** on line 6 with value 1.40. Only one item (the **K** on line 6) was missed prior to that, so only  $1 \times 0.04$  is subtracted as the scoring correction. The final scored log contrast sensitivity is thus 1.36.

Row	FORM 1	Left eye <input checked="" type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>				
1	C <input type="checkbox"/> 0.04	H <input type="checkbox"/> 0.08	V <input type="checkbox"/> 0.12	O <input type="checkbox"/> 0.16	S <input type="checkbox"/> 0.20	N <input type="checkbox"/> 0.24	Enter value at final correct item	1.40
2	D <input type="checkbox"/> 0.28	S <input type="checkbox"/> 0.32	Z <input type="checkbox"/> 0.36	N <input type="checkbox"/> 0.40	R <input type="checkbox"/> 0.44	K <input type="checkbox"/> 0.48	Enter number of errors prior to final correct item	2 <u>1</u>
3	N <input type="checkbox"/> 0.52	D <input type="checkbox"/> 0.56	R <input type="checkbox"/> 0.60	H <input type="checkbox"/> 0.64	V <input type="checkbox"/> 0.68	Z <input type="checkbox"/> 0.72	Multiply line 2 by 0.04	3 <u>0.04</u>
4	C <input type="checkbox"/> 0.76	S <input type="checkbox"/> 0.80	O <input type="checkbox"/> 0.84	N <input type="checkbox"/> 0.88	K <input type="checkbox"/> 0.92	H <input type="checkbox"/> 0.96		
5	K <input type="checkbox"/> 1.00	N <input type="checkbox"/> 1.04	V <input type="checkbox"/> 1.08	D <input type="checkbox"/> 1.12	S <input type="checkbox"/> 1.16	R <input type="checkbox"/> 1.20		
6	Z <input type="checkbox"/> 1.24	R <input type="checkbox"/> 1.28	D <input type="checkbox"/> 1.32	K <input checked="" type="checkbox"/> 1.36	H <input type="checkbox"/> 1.40	O <input checked="" type="checkbox"/> 1.44	Subtract line 3 from line 1	
7	H <input checked="" type="checkbox"/> 1.48	Z <input type="checkbox"/> 1.52	C <input type="checkbox"/> 1.56	V <input type="checkbox"/> 1.60	R <input type="checkbox"/> 1.64	K <input type="checkbox"/> 1.68	Log CS Score:	1.36
8	S <input type="checkbox"/> 1.72	C <input type="checkbox"/> 1.76	Z <input type="checkbox"/> 1.80	D <input type="checkbox"/> 1.84	V <input type="checkbox"/> 1.88	O <input type="checkbox"/> 1.92		

## Score Interpretation

The table below may be used as a rough guide for interpreting log CS scores in relation to severity of contrast sensitivity loss, with the middle column showing severity ranges in terms of test score, and the right column showing severity in terms of which lines the examinee can and cannot read.

Severity	log CS score range	Line-based range
<b>Profound loss</b>	$\log CS < 0.48$	Cannot read to the end of line 2
<b>Severe loss</b>	$0.52 \leq \log CS \leq 1.00$	Reads lines 1–2, but cannot read all of lines 3 and 4
<b>Moderate loss</b>	$1.04 \leq \log CS \leq 1.48$	Reads lines 1–4 but cannot read beyond line 6
<b>Mild loss</b> (youth and adults $\leq$ age 60) <b>No loss</b> ( $>$ age 60)	$1.52 \leq \log CS \leq 1.76$	Reads lines 1–6 but cannot read beyond beginning of line 7
<b>No loss</b>	$1.72 \leq \log CS \leq 1.92$	Reads beyond lines 1-7

Note that when the two eyes have similar log CS scores, binocular log CS scores may be expected to be about 0.15 higher than monocular log CS scores.




Study.<sup>2</sup> Both tests have optotypes arranged in eight rows of six letters. Visual size of the PRCSC letters at 1 m test distance is 2.8 deg, similar to the 2.5 deg letters of the Mars at 40 cm. Both tests also assess true contrast sensitivity.

The Mars test, however, improves on the PRCSC in six important ways. (1) The Mars test is small, measuring 24 x 36 cm (9.5 x 14 inches), and can be held in the hands, whereas the PRCSC measures 59 cm x 84 cm (23 x 33 inches), and requires wall mounting. (2) The Mars test has 0.04 log unit contrast differences from letter-to-letter (or numeral-to-numeral), as compared to the PRCSC, which has 0.15 log unit gradations. Since the Mars test's scores take on values to the nearest 0.04 log unit, finer differences in contrast sensitivity may be discriminated. (3) The Mars test can be stored and protected after use, whereas the PRCSC is vulnerable to bruises, nicks, dust, and deterioration from constant exposure to light. (4) The Mars test requires even illumination of only 858 cm<sup>2</sup>, compared to the 4,956 cm<sup>2</sup> that must be illuminated on the PRCSC (9,912 cm<sup>2</sup> if both PRCSC charts are mounted). (5) The recommended Mars test distance is the same as the one used for standard clinical refraction for reading, so the test can be administered while the patient is already in the examination chair and properly corrected. The PRCSC requires re-seating the examinee to a special distance of 1 m (3 ft) distance (from the wall holding the chart). This is a distance that presbyopic patients are not generally corrected to, and so without the use of additional lenses, the chart image may be slightly out of focus during testing. (6) The Mars test is portable and can be used for testing in the home or other remote testing sites.

Because the PRCSC was the first letter contrast sensitivity test, and was used in early studies, many researchers assume that using it will facilitate comparison with other studies. But over the years, the PRCSC has undergone many changes that invalidate this assumption. For example, the recommended test distance is now 1 m (3 ft), one-third of the 3 m (10 ft) distance originally recommended.<sup>5</sup> The PRCSC also now tests a sensitivity range that is shifted by 0.05 log unit relative to the original test.<sup>5</sup> But most important are the many changes to test procedures adopted by researchers and other users seeking to increase the test's accuracy. The PRCSC presents letter optotypes in triplets of the same contrast. The manufacturer's instructions state that testing stops when two or more errors are made within a triplet, and the final score is the log CS value of the lowest contrast triplet in which at least two of the letters has been correctly identified. But many (if not most) researchers substitute so called 'letter-by-letter' scoring that gives credit for correct identifications within a triplet.<sup>7</sup> Additionally, many alter the stopping rule. Some stop when there are two or more errors, while others stop the test when all three identifications within a triplet are incorrect. A further complication is that some researchers recommend allowing a response of "O" for a presented "C" to be scored as correct,<sup>8</sup> convoluting interpretation of the guessing probability. These may seem like minor details, but guessing probability, and rules for stopping and scoring are all known to have a big impact on test scores.<sup>4,5</sup> The important point is that the wide variation in PRCSC administration and scoring make it difficult to interpret population data or develop consistent population norms.<sup>4</sup> Because of all these variations in usage, there is no single PRCSC test that has been adopted by consensus in the research community. The Mars test, on the other hand is a single test with a single scoring method, that is inherently more accurate than the PRCSC, as shown in the following analysis.

One of the most important parameters that determines design accuracy of tests like the Mars and PRCSC is the so-called misreporting rate,<sup>4,5</sup> also known as the lapse rate. This parameter is the frequency at which those being tested will misspeak and report a letter other than the one they saw. The original designers of the PRCSC made the assumption (perhaps thinking that verbal reporting would be similar to computer-controlled testing where subjects respond by pressing buttons) that people will misreport at a rate ranging from 1% to 5%. Measurement of the rate using the verbal letter identification task used in letter contrast sensitivity tests however, has shown the PRCSC's design assumption to be off by more than an order of magnitude. The empirically

estimated rate is actually only about 0.05%.<sup>9</sup> The graph at the right shows how test accuracy decreases with higher misreporting rates for both the Mars and the PRCSC (note that lower accuracy is associated with *higher* values on the vertical axis). Both tests show accuracy start to decline rapidly when misreporting is just below 1%. It also shows that both tests have similar, but relatively poor accuracy in the 1–5% misreporting range. But most importantly, it shows that at the negligible rates of misreporting empirically observed using the verbal responses elicited in chart-based testing, the Mars tests (with their much finer contrast gradations and superior scoring method) can be expected to be 28% more accurate than the PRCSC.



Thus, the Mars test has improved accuracy and simpler scoring, as well as being more convenient to use in clinical settings.

## Specifications

Each Mars chart measures 9.5 x 14 inches (24.1 x 35.6 cm). Optotypes are composed within a 1.75 cm square box, which at 15.75 inches (40 cm) viewing distance, subtends 2.5 deg of visual angle. Each form consists of 48 optotypes, arranged in eight rows of six. The contrast of each optotype, reading from left to right, and continuing on successive lines on each chart, decreases by a constant factor of 0.04 log unit, or 0.912 from the reference contrast of 1.0. The highest contrast optotype on each of the Mars chart forms (at row 1, column 1) has a contrast of 0.912, and the lowest (at row 8, column 6), has a contrast of 0.012. The table below shows the contrast values of all the letters or numerals on each Mars form. Contrast ( $C$ ) is defined as

$$C = \frac{R_{fg} - R_{bg}}{R_{bg}},$$

where  $R_{fg}$  and  $R_{bg}$  are the reflectances of the image foreground (the printed letter or numeral) and background (the white background), respectively. If threshold contrast is  $C$ , then log contrast sensitivity ( $\log CS$ ) associated with each contrast value is  $-\log_{10} C$ .

Column⇒	1		2		3		4		5		6	
Row ↓	$\log CS$	$C$										
1	0.04	0.912	0.08	0.832	0.12	0.759	0.16	0.692	0.20	0.631	0.24	0.575
2	0.28	0.525	0.32	0.479	0.36	0.437	0.40	0.398	0.44	0.363	0.48	0.331
3	0.52	0.302	0.56	0.275	0.60	0.251	0.64	0.229	0.68	0.209	0.72	0.191
4	0.76	0.174	0.80	0.158	0.84	0.145	0.88	0.132	0.92	0.120	0.96	0.110
5	1.00	0.100	1.04	0.091	1.08	0.083	1.12	0.076	1.16	0.069	1.20	0.063
6	1.24	0.058	1.28	0.052	1.32	0.048	1.36	0.044	1.40	0.040	1.44	0.036
7	1.48	0.033	1.52	0.030	1.56	0.028	1.60	0.025	1.64	0.023	1.68	0.021
8	1.72	0.019	1.76	0.017	1.80	0.016	1.84	0.014	1.88	0.013	1.92	0.012

All three forms in each version of the test have the same set of contrasts, with only the letters or numerals on each chart varying. For example, the contrast of the optotype in row 4, column 2 is 0.158. Reading correctly only to that letter is associated with a log CS score of 0.80.

## How to Cite the Mars Test in Published Works

If you use the Mars Test in your research or other publications, please cite Reference 4 below (shown in **bold**).

### Literature Cited

1. National Research Council (U.S.), Committee on Vision, Working Group 39. Recommended standard procedures for the clinical measurement and specification of visual acuity. Basel; New York: S. Karger; 1980.
2. Ferris FL, Kassoff A, Bresnick GH, Bailey I. New Visual Acuity Charts for Clinical Research. *Am J Ophthalmol* [Internet]. 1982 Jul 1;94(1):91–96. Available from: <http://www.sciencedirect.com/science/article/pii/0002939482901970>
3. Rosenbaum R, Fleischmann M. Character Frequency in Multilingual Corpus 1 – Part 1. *J Quant Linguist* [Internet]. 2002 Dec 1;9(3):233–260. Available from: <https://doi.org/10.1076/jqul.9.3.233.14122>
4. **Arditi A. Improving the design of the letter contrast sensitivity test. *Invest Ophthalmol Vis Sci* [Internet]. 2005;46(6):2225–2229. Available from: <https://iovs.arvojournals.org/article.aspx?articleid=2163640>**
5. Pelli D, Robson J, Wilkins, AJ. The design of a new letter chart for measuring contrast sensitivity. *Clin Vis Sci*. 1988;2:187–199.
6. Sloan LL. New test Charts for the Measurement of Visual Acuity at far and Near Distances. *Am J Ophthalmol* [Internet]. 1959 Dec 1 [cited 2018 Nov 17];48(6):807–813. Available from: [https://doi.org/10.1016/0002-9394\(59\)90626-9](https://doi.org/10.1016/0002-9394(59)90626-9)
7. Elliott DB, Sanderson K, Conkey A. The reliability of the Pelli-Robson contrast sensitivity chart. *Ophthalmic Physiol Opt*. 1990;10(1):21–24.
8. Elliott DB, Whitaker D, Bonette L. Differences in the legibility of letters at contrast threshold using the Pelli-Robson chart. *Ophthalmic Physiol Opt* [Internet]. 1990 Oct 1 [cited 2018 Dec 15];10(4):323–326. Available from: <https://doi.org/10.1111/j.1475-1313.1990.tb00877.x>
9. Arditi A. Lapse resistance in the verbal letter reporting task. *Vision Res*. 2006;46(8–9):1327–1330.

# The Mars Letter Contrast Sensitivity Test

## Score Sheet

Examinee \_\_\_\_\_ Date \_\_\_\_\_

Correction \_\_\_\_\_ Test distance \_\_\_\_\_ Administered by \_\_\_\_\_

Comments \_\_\_\_\_

Instruct examinee to read letters left to right for each line, from top to bottom of the chart

**Important: Allow *only* the letters C D H K N O R S V Z as responses.**

Mark errors with an "X." Stop testing after **two** consecutive errors.

Row	FORM 1	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	C <input type="checkbox"/> 0.04	H <input type="checkbox"/> 0.08	V <input type="checkbox"/> 0.12	O <input type="checkbox"/> 0.16	S <input type="checkbox"/> 0.20 N <input type="checkbox"/> 0.24	Enter value at final correct item 1 _____
2	D <input type="checkbox"/> 0.28	S <input type="checkbox"/> 0.32	Z <input type="checkbox"/> 0.36	N <input type="checkbox"/> 0.40	R <input type="checkbox"/> 0.44 K <input type="checkbox"/> 0.48	Enter number of errors prior to final correct item 2 _____
3	N <input type="checkbox"/> 0.52	D <input type="checkbox"/> 0.56	R <input type="checkbox"/> 0.60	H <input type="checkbox"/> 0.64	V <input type="checkbox"/> 0.68 Z <input type="checkbox"/> 0.72	Multiply line 2 by 0.04 3 _____
4	C <input type="checkbox"/> 0.76	S <input type="checkbox"/> 0.80	O <input type="checkbox"/> 0.84	N <input type="checkbox"/> 0.88	K <input type="checkbox"/> 0.92 H <input type="checkbox"/> 0.96	Subtract line 3 from line 1 <b>Log CS Score: _____</b>
5	K <input type="checkbox"/> 1.00	N <input type="checkbox"/> 1.04	V <input type="checkbox"/> 1.08	D <input type="checkbox"/> 1.12	S <input type="checkbox"/> 1.16 R <input type="checkbox"/> 1.20	
6	Z <input type="checkbox"/> 1.24	R <input type="checkbox"/> 1.28	D <input type="checkbox"/> 1.32	K <input type="checkbox"/> 1.36	H <input type="checkbox"/> 1.40 O <input type="checkbox"/> 1.44	
7	H <input type="checkbox"/> 1.48	Z <input type="checkbox"/> 1.52	C <input type="checkbox"/> 1.56	V <input type="checkbox"/> 1.60	R <input type="checkbox"/> 1.64 K <input type="checkbox"/> 1.68	
8	S <input type="checkbox"/> 1.72	C <input type="checkbox"/> 1.76	Z <input type="checkbox"/> 1.80	D <input type="checkbox"/> 1.84	V <input type="checkbox"/> 1.88 O <input type="checkbox"/> 1.92	

Row	FORM 2	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	K <input type="checkbox"/> 0.04	S <input type="checkbox"/> 0.08	H <input type="checkbox"/> 0.12	O <input type="checkbox"/> 0.16	N <input type="checkbox"/> 0.20 C <input type="checkbox"/> 0.24	Enter value at final correct item 1 _____
2	Z <input type="checkbox"/> 0.28	D <input type="checkbox"/> 0.32	C <input type="checkbox"/> 0.36	R <input type="checkbox"/> 0.40	V <input type="checkbox"/> 0.44 O <input type="checkbox"/> 0.48	Enter number of errors prior to final correct item 2 _____
3	C <input type="checkbox"/> 0.52	K <input type="checkbox"/> 0.56	O <input type="checkbox"/> 0.60	N <input type="checkbox"/> 0.64	R <input type="checkbox"/> 0.68 S <input type="checkbox"/> 0.72	Multiply line 2 by 0.04 3 _____
4	N <input type="checkbox"/> 0.76	S <input type="checkbox"/> 0.80	Z <input type="checkbox"/> 0.84	K <input type="checkbox"/> 0.88	H <input type="checkbox"/> 0.92 D <input type="checkbox"/> 0.96	Subtract line 3 from line 1 <b>Log CS Score: _____</b>
5	H <input type="checkbox"/> 1.00	N <input type="checkbox"/> 1.04	C <input type="checkbox"/> 1.08	O <input type="checkbox"/> 1.12	R <input type="checkbox"/> 1.16 Z <input type="checkbox"/> 1.20	
6	V <input type="checkbox"/> 1.24	K <input type="checkbox"/> 1.28	S <input type="checkbox"/> 1.32	N <input type="checkbox"/> 1.36	D <input type="checkbox"/> 1.40 R <input type="checkbox"/> 1.44	
7	K <input type="checkbox"/> 1.48	R <input type="checkbox"/> 1.52	V <input type="checkbox"/> 1.56	Z <input type="checkbox"/> 1.60	O <input type="checkbox"/> 1.64 S <input type="checkbox"/> 1.68	
8	V <input type="checkbox"/> 1.72	Z <input type="checkbox"/> 1.76	C <input type="checkbox"/> 1.80	D <input type="checkbox"/> 1.84	V <input type="checkbox"/> 1.88 H <input type="checkbox"/> 1.92	

Row	FORM 3	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	H <input type="checkbox"/> 0.04	R <input type="checkbox"/> 0.08	Z <input type="checkbox"/> 0.12	V <input type="checkbox"/> 0.16	C <input type="checkbox"/> 0.20 N <input type="checkbox"/> 0.24	Enter value at final correct item 1 _____
2	S <input type="checkbox"/> 0.28	O <input type="checkbox"/> 0.32	K <input type="checkbox"/> 0.36	D <input type="checkbox"/> 0.40	R <input type="checkbox"/> 0.44 S <input type="checkbox"/> 0.48	Enter number of errors prior to final correct item 2 _____
3	K <input type="checkbox"/> 0.52	D <input type="checkbox"/> 0.56	C <input type="checkbox"/> 0.60	V <input type="checkbox"/> 0.64	O <input type="checkbox"/> 0.68 H <input type="checkbox"/> 0.72	Multiply line 2 by 0.04 3 _____
4	N <input type="checkbox"/> 0.76	S <input type="checkbox"/> 0.80	O <input type="checkbox"/> 0.84	Z <input type="checkbox"/> 0.88	C <input type="checkbox"/> 0.92 D <input type="checkbox"/> 0.96	Subtract line 3 from line 1 <b>Log CS Score: _____</b>
5	R <input type="checkbox"/> 1.00	H <input type="checkbox"/> 1.04	N <input type="checkbox"/> 1.08	K <input type="checkbox"/> 1.12	Z <input type="checkbox"/> 1.16 O <input type="checkbox"/> 1.20	
6	C <input type="checkbox"/> 1.24	R <input type="checkbox"/> 1.28	S <input type="checkbox"/> 1.32	V <input type="checkbox"/> 1.36	K <input type="checkbox"/> 1.40 N <input type="checkbox"/> 1.44	
7	S <input type="checkbox"/> 1.48	K <input type="checkbox"/> 1.52	R <input type="checkbox"/> 1.56	N <input type="checkbox"/> 1.60	H <input type="checkbox"/> 1.64 D <input type="checkbox"/> 1.68	
8	C <input type="checkbox"/> 1.72	V <input type="checkbox"/> 1.76	H <input type="checkbox"/> 1.80	D <input type="checkbox"/> 1.84	O <input type="checkbox"/> 1.88 Z <input type="checkbox"/> 1.92	

# The Mars Numeral Contrast Sensitivity Test

## Score Sheet

Examinee \_\_\_\_\_ Date \_\_\_\_\_

Correction \_\_\_\_\_ Test distance \_\_\_\_\_ Administered by \_\_\_\_\_

Comments \_\_\_\_\_

Instruct examinee to read numerals left to right for each line, from top to bottom of the chart

**Important: Allow *only* the numerals **0 1 2 3 4 5 6 7 8 9** as responses.**

Mark errors with an "X." Stop testing after **two** consecutive errors.

Row	FORM 1	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	0 <input type="checkbox"/> 0.04	2 <input type="checkbox"/> 0.08	8 <input type="checkbox"/> 0.12	5 <input type="checkbox"/> 0.16	7 <input type="checkbox"/> 0.20	4 <input type="checkbox"/> 0.24
2	1 <input type="checkbox"/> 0.28	7 <input type="checkbox"/> 0.32	9 <input type="checkbox"/> 0.36	4 <input type="checkbox"/> 0.40	6 <input type="checkbox"/> 0.44	3 <input type="checkbox"/> 0.48
3	4 <input type="checkbox"/> 0.52	1 <input type="checkbox"/> 0.56	6 <input type="checkbox"/> 0.60	2 <input type="checkbox"/> 0.64	8 <input type="checkbox"/> 0.68	9 <input type="checkbox"/> 0.72
4	0 <input type="checkbox"/> 0.76	7 <input type="checkbox"/> 0.80	5 <input type="checkbox"/> 0.84	4 <input type="checkbox"/> 0.88	3 <input type="checkbox"/> 0.92	2 <input type="checkbox"/> 0.96
5	3 <input type="checkbox"/> 1.00	4 <input type="checkbox"/> 1.04	8 <input type="checkbox"/> 1.08	1 <input type="checkbox"/> 1.12	7 <input type="checkbox"/> 1.16	6 <input type="checkbox"/> 1.20
6	9 <input type="checkbox"/> 1.24	6 <input type="checkbox"/> 1.28	1 <input type="checkbox"/> 1.32	3 <input type="checkbox"/> 1.36	2 <input type="checkbox"/> 1.40	5 <input type="checkbox"/> 1.44
7	2 <input type="checkbox"/> 1.48	9 <input type="checkbox"/> 1.52	0 <input type="checkbox"/> 1.56	8 <input type="checkbox"/> 1.60	6 <input type="checkbox"/> 1.64	3 <input type="checkbox"/> 1.68
8	7 <input type="checkbox"/> 1.72	0 <input type="checkbox"/> 1.76	9 <input type="checkbox"/> 1.80	1 <input type="checkbox"/> 1.84	8 <input type="checkbox"/> 1.88	5 <input type="checkbox"/> 1.92

Enter value at final correct item 1 \_\_\_\_\_

Enter number of errors prior to final correct item 2 \_\_\_\_\_

Multiply line 2 by 0.04 3 \_\_\_\_\_

Subtract line 3 from line 1  
Log CS Score: \_\_\_\_\_

Row	FORM 2	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	3 <input type="checkbox"/> 0.04	7 <input type="checkbox"/> 0.08	2 <input type="checkbox"/> 0.12	5 <input type="checkbox"/> 0.16	4 <input type="checkbox"/> 0.20	0 <input type="checkbox"/> 0.24
2	9 <input type="checkbox"/> 0.28	1 <input type="checkbox"/> 0.32	0 <input type="checkbox"/> 0.36	6 <input type="checkbox"/> 0.40	8 <input type="checkbox"/> 0.44	5 <input type="checkbox"/> 0.48
3	0 <input type="checkbox"/> 0.52	3 <input type="checkbox"/> 0.56	5 <input type="checkbox"/> 0.60	4 <input type="checkbox"/> 0.64	6 <input type="checkbox"/> 0.68	7 <input type="checkbox"/> 0.72
4	4 <input type="checkbox"/> 0.76	7 <input type="checkbox"/> 0.80	9 <input type="checkbox"/> 0.84	3 <input type="checkbox"/> 0.88	2 <input type="checkbox"/> 0.92	1 <input type="checkbox"/> 0.96
5	2 <input type="checkbox"/> 1.00	4 <input type="checkbox"/> 1.04	0 <input type="checkbox"/> 1.08	5 <input type="checkbox"/> 1.12	6 <input type="checkbox"/> 1.16	9 <input type="checkbox"/> 1.20
6	8 <input type="checkbox"/> 1.24	3 <input type="checkbox"/> 1.28	7 <input type="checkbox"/> 1.32	4 <input type="checkbox"/> 1.36	1 <input type="checkbox"/> 1.40	6 <input type="checkbox"/> 1.44
7	3 <input type="checkbox"/> 1.48	6 <input type="checkbox"/> 1.52	8 <input type="checkbox"/> 1.56	9 <input type="checkbox"/> 1.60	5 <input type="checkbox"/> 1.64	7 <input type="checkbox"/> 1.68
8	8 <input type="checkbox"/> 1.72	9 <input type="checkbox"/> 1.76	0 <input type="checkbox"/> 1.80	1 <input type="checkbox"/> 1.84	8 <input type="checkbox"/> 1.88	2 <input type="checkbox"/> 1.92

Enter value at final correct item 1 \_\_\_\_\_

Enter number of errors prior to final correct item 2 \_\_\_\_\_

Multiply line 2 by 0.04 3 \_\_\_\_\_

Subtract line 3 from line 1  
Log CS Score: \_\_\_\_\_

Row	FORM 3	Left eye <input type="checkbox"/>	Right eye <input type="checkbox"/>	Binocular <input type="checkbox"/>		
1	2 <input type="checkbox"/> 0.04	6 <input type="checkbox"/> 0.08	9 <input type="checkbox"/> 0.12	8 <input type="checkbox"/> 0.16	0 <input type="checkbox"/> 0.20	4 <input type="checkbox"/> 0.24
2	7 <input type="checkbox"/> 0.28	5 <input type="checkbox"/> 0.32	3 <input type="checkbox"/> 0.36	1 <input type="checkbox"/> 0.40	6 <input type="checkbox"/> 0.44	7 <input type="checkbox"/> 0.48
3	3 <input type="checkbox"/> 0.52	1 <input type="checkbox"/> 0.56	0 <input type="checkbox"/> 0.60	8 <input type="checkbox"/> 0.64	5 <input type="checkbox"/> 0.68	2 <input type="checkbox"/> 0.72
4	4 <input type="checkbox"/> 0.76	7 <input type="checkbox"/> 0.80	5 <input type="checkbox"/> 0.84	9 <input type="checkbox"/> 0.88	0 <input type="checkbox"/> 0.92	1 <input type="checkbox"/> 0.96
5	6 <input type="checkbox"/> 1.00	2 <input type="checkbox"/> 1.04	4 <input type="checkbox"/> 1.08	3 <input type="checkbox"/> 1.12	9 <input type="checkbox"/> 1.16	5 <input type="checkbox"/> 1.20
6	0 <input type="checkbox"/> 1.24	6 <input type="checkbox"/> 1.28	7 <input type="checkbox"/> 1.32	8 <input type="checkbox"/> 1.36	3 <input type="checkbox"/> 1.40	4 <input type="checkbox"/> 1.44
7	7 <input type="checkbox"/> 1.48	3 <input type="checkbox"/> 1.52	6 <input type="checkbox"/> 1.56	4 <input type="checkbox"/> 1.60	2 <input type="checkbox"/> 1.64	1 <input type="checkbox"/> 1.68
8	0 <input type="checkbox"/> 1.72	8 <input type="checkbox"/> 1.76	2 <input type="checkbox"/> 1.80	1 <input type="checkbox"/> 1.84	5 <input type="checkbox"/> 1.88	9 <input type="checkbox"/> 1.92

Enter value at final correct item 1 \_\_\_\_\_

Enter number of errors prior to final correct item 2 \_\_\_\_\_

Multiply line 2 by 0.04 3 \_\_\_\_\_

Subtract line 3 from line 1  
Log CS Score: \_\_\_\_\_