

The GI Effects Comprehensive Report Review: A Global Evaluation of GI Function

Stephen L. Goldman, DC
Genova Diagnostics





Stephen L. Goldman, DC

Medical Education Specialist



GEN
DIAG

Patient:

2200 GI

Key

Dige

MAL

Therapeutic Support Options

Product Breakdown
Fecal F
Pancrea
• Digest
• Betain
• Bile Sa
• Apple
• Mindu
• Digest

Commensal

Patient Total

Total Commensal
healthy cohort.
prebiotic-rich for
potential bacter


Dysbiosis

Inflammation

Methane Dys

Immune Suppression
(Metabolic Score)

Commensal



Relative

Bacteroid

Firmicute

Actinoba

Proteoba

Euryarch

Fusobac

Verrucom

Relative

indicate b

appropria

nutrient in

Approxim

Physi

2200 GI Effects™

Methodology: GC-FID, Auto

Pancreatic Elastase 1 †

Products of Protein Breakdown
(Valerate, Isobutyrate)

Fecal Fat (Total*)

Triglycerides

Long-Chain Fatty Acids

Cholesterol

Phospholipids

Calprotectin †

Eosinophil Protein X (EPX)

Fecal secretory IgA

Metabolic

Short-Chain Fatty Acids
(Acetate, n-Butyrate)

n-Butyrate Concentration

n-Butyrate %

Acetate %

Propionate %

Beta-glucuronidase

*Total value is equal to the sum of the individual values.

†These results are not reported.

Tests were developed and validated for use with the 2200 GI Effects™.

Methodology: DNA by qPCR

Commensal Bacteria

Bacteroidetes Phylum

Bacteroides uniformis

Phocaeicola vulgatus

Barnesiella spp.

Odonibacter spp.

Prevotella spp.

Firmicutes Phylum

Anaerotruncus colih

Butyrivibrio crossotus

Clostridium spp.

Coprococcus eutacti

Faecalibacterium pr

Lactobacillus spp.

Pseudoflavonifractor

Roseburia spp.

Ruminococcus brom

Veillonella spp.

Actinobacteria Phylum

Bifidobacterium spp.

Bifidobacterium

Collinsella aerofacie

Proteobacteria Phylum

Desulfovibrio piger

Escherichia coli

Oxalobacter formige

Euryarchaeota Phylum

Methanobrevibacter

Fusobacteria Phylum

Fusobacterium spp.

Methodology: Culture/MAL

Human microflora is a competitive ecosystem. The significance should be determined by the reference of all patients.

NG

No Growth

Bacteriology

Lactobacillus spp.

Escherichia coli

Bifidobacterium (Anaerobes)

Additional Bacteria

Klebsiella pneumoniae

Bacillus species

Enterococcus faecium

Mycology (Culture)

Candida krusei

Yeast, not Candida albicans

Methodology: Potassium Hydroxide

Potassium Hydroxide

These yeast usually grow on potassium hydroxide. However, moderate to strong growth is not observed.

KOH Preparation

One negative specimen

Methodology: Vitek

Microscopic O&P

Microscopic O&P is commonly found in the Additional Findings section of all patients.

Genus/species

Nematodes - roundworms

Ancylostoma/Necator

Ascaris lumbricoides

Capillaria philippinensis

Enterobius vermiciformis

Strongyloides stercoralis

Trichuris trichiura

Cestodes - tapeworms

Diphyllobothrium latum

Dipylidium caninum

Hymenolepis diminuta

Hymenolepis nana

Taenia spp.

Trematodes - flukes

Clonorchis/Opisthorchis

Fasciola spp./ Fasciolopsis

Heterophyes/Metacercariae

Paragonimus spp.

Schistosoma spp.

Protozoa

Balantidium coli

Blastocystis spp.

Chilomastix mesnili

Cryptosporidium spp.

Cyclospora cayentensis

Dientamoeba fragilis

Entamoeba coli

Methodology: EIA

PCR Parasitology

Organism

Blastocystis spp.

Cryptosporidium parvum

Cyclospora cayentensis

Dientamoeba fragilis

Entamoeba histolytica

Giardia

Methodology: Fecal Immunochemical Test

Fecal Occult Blood*

Color††

Consistency††

††Results provided from the 2200 GI Effects™. Tests were developed and validated for use with the 2200 GI Effects™.

Methodology: EIA

Zonulin Family Peptide

Reference:

1. Scheffler L, et al. WJG. 2013;27(12):1381-1386.

Recognizes Properdin

Methodology: Mac

No human parasite

Methodology: Vitek

Klebsiella pneumoniae

Ampicillin

Amox/Clav

Cephalexin

Ciprofloxacin

Tetracycline

Trimethoprim

Natural Agents

Klebsiella pneumoniae

Berberine

Oregano

Uva-Ursi

Methodology: Vitek

HpSA - H. pylori

Campylobacter

Clostridium difficile

Shiga toxin E.

Fecal Lactoferrin

Prescriptive Agents

The R (Resistant)

The I (Intermediate)

levels and for which the S-DD (Susceptible) was not achieved.

The S (Susceptible)

NI (No Interpretation)

Refer to published literature for further information.

Natural Agents

In this assay, inhibition is defined as the reduction level on organism growth as a direct result of inhibition by a natural substance. The level of inhibition is an indicator of how effective the substance was at limiting the growth of an organism in an in vitro environment. High inhibition indicates a greater ability by the substance to limit growth, while Low Inhibition a lesser ability to limit growth. The designated natural products should be considered investigational in nature and not be viewed as standard clinical treatment substances.

Tests were developed and validated for use with the 2200 GI Effects™.

Methodology: Vitek

Prescriptive Agents

Klebsiella pneumoniae

Ampicillin

Amox/Clav

Cephalexin

Ciprofloxacin

Tetracycline

Trimethoprim

Natural Agents

Klebsiella pneumoniae

Berberine

Oregano

Uva-Ursi

Methodology: Vitek

HpSA - H. pylori

Campylobacter

Clostridium difficile

Shiga toxin E.

Fecal Lactoferrin

Prescriptive Agents

The R (Resistant)

The I (Intermediate)

levels and for which the S-DD (Susceptible) was not achieved.

The S (Susceptible)

NI (No Interpretation)

Refer to published literature for further information.

Natural Agents

In this assay, inhibition is defined as the reduction level on organism growth as a direct result of inhibition by a natural substance. The level of inhibition is an indicator of how effective the substance was at limiting the growth of an organism in an in vitro environment. High inhibition indicates a greater ability by the substance to limit growth, while Low Inhibition a lesser ability to limit growth. The designated natural products should be considered investigational in nature and not be viewed as standard clinical treatment substances.

Tests were developed and validated for use with the 2200 GI Effects™.

Methodology: Vitek 2® System Microbial Antibiotic susceptibility, Manual Minimum Inhibition Concentration

Mycology Sensitivity

Candida Susceptibility Profile for Azoles*

Organism	Number of Isolates	% Sensitive	
		Fluconazole	Voriconazole
Candida albicans	25561	99.19%	99.51%
Candida parapsilosis	8777	98.64%	99.33%
Candida krusei	3420	0.23%	97.79%
Candida tropicalis	1076	93.22%	90.57%
Candida glabrata	2898	27.1%	90.9%

*Results of pharmaceutical sensitivities against certain yeast species are based on internal Genova data pertaining to the frequency of susceptibility of the specific yeast to the listed antifungal agent. The pharmaceutical results are not patient-specific. Conversely, the results of inhibition to nystatin and natural agents are patient-specific.

Non-absorbed Antifungals

Candida krusei	LOW INHIBITION	HIGH INHIBITION
Nystatin		

Natural Agents

Candida krusei	LOW INHIBITION	HIGH INHIBITION
Berberine		
Caprylic Acid		
Garlic		
Undecylenic Acid		
Uva-Ursi		

Nystatin and Natural Agents:

Results for Nystatin are being reported with natural antifungals in this category in accordance with laboratory guidelines for reporting sensitivities. In this assay, inhibition is defined as the reduction level on organism growth as a direct result of inhibition by a natural substance. The level of inhibition is an indicator of how effective the substance was at limiting the growth of an organism in an in vitro environment. High inhibition indicates a greater ability by the substance to limit growth, while Low Inhibition a lesser ability to limit growth. The designated natural products should be considered investigational in nature and not be viewed as standard clinical treatment substances.

Therapeutic Support



GENOVA
DIAGNOSTICS

63 Zillicoa Street
Asheville, NC 28801
© Genova Diagnostics



Patient:

2200 GI Effects™ Comprehensive Profile - Stool

Powered by **Genova AI**

Results Overview




Functional Imbalance Scores

Key (<2): Low Need for Support (2-3): Optional Need for Support (4-6): Moderate Need for Support (7-10): High Need for Support

	Need for Digestive Support	Need for Inflammation Modulation	Need for Microbiome Support	Need for Prebiotic Support	Need for Antimicrobial Support
	MALDIGESTION	INFLAMMATION	DYSBIOSIS	METABOLIC IMBALANCE	INFECTION
	0	0	4	2	7
Biomarkers	<ul style="list-style-type: none"> Products of Protein Breakdown Fecal Fats Pancreatic Elastase 	<ul style="list-style-type: none"> Secretory IgA Calprotectin Eosinophil Protein X Occult Blood 	<ul style="list-style-type: none"> PP Bacteria/Yeast IAD/Methane Score Reference Variance Total Abundance 	<ul style="list-style-type: none"> Total SCFA's n-Butyrate Conc. SCFA (%) Beta-glucuronidase 	<ul style="list-style-type: none"> PP Bacteria/Yeast Parasitic Infection Pathogenic Bacteria Total Abundance
Therapeutic Support Options	<ul style="list-style-type: none"> Digestive Enzymes Betaine HCl Bile Salts Apple Cider Vinegar Mindful Eating Habits Digestive Bitters 	<ul style="list-style-type: none"> Elimination Diet/ Food Sensitivity Testing Mucosa Support: Slippery Elm, Althea, Aloe, DGL, etc. Zinc Carnosine L-Glutamine Quercetin Turmeric Omega-3's GI Referral (if Calpro is Elevated) 	<ul style="list-style-type: none"> Pre-/Probiotics Increase Dietary Fiber Intake Consider SIBO Testing Increase Resistant Starches Increase Fermented Foods Meal Timing 	<ul style="list-style-type: none"> Pre-/Probiotics Increase Dietary Fiber Intake Increase Resistant Starches Increase Fermented Foods Calcium D-Glucarate (for high beta-glucuronidase) 	<ul style="list-style-type: none"> Antibiotics (if warranted) Antimicrobial Herbal Therapy Antiparasitic Herbal Therapy (if warranted) Saccharomyces boulardii

ons left up to clinician

Commensal Bacteria




- Clinicians often struggle with what to do with DNA PCR analysis of commensal bacteria
- Historical limitations
 - Methodologies differ in literature
 - Discrepant results in publications
 - Unknown clinical importance of individual bacteria
 - No research into bacterial patterns

A Novel Approach to Microbiome Analysis


1. Abundance

2. Patterns

3. Balance



Total Commensal Balance: The total commensal abundance is a sum-total of the reported commensal bacteria compared to a healthy cohort. Low levels of commensal bacteria are often observed after antimicrobial therapy, or in diets lacking fiber and/or prebiotic-rich foods and may indicate the need for microbiome support. Conversely, higher total commensal abundance may indicate






Dysbiosis Patterns: Genova's data analysis has led to the development of unique dysbiosis patterns, related to key physiologic disruptions, such as immunosuppression and inflammation. These patterns may represent dysbiotic changes that could pose clinical significance. Please see Genova's published literature for more details: <https://dcu.be/bRhzv>

Zone 1: The commensal profile in this zone does not align with profiles associated with intestinal inflammation or immunosuppression. If inflammatory biomarkers are present, other causes need to be excluded, such as infection, food allergy, or more serious pathology.

Zone 2: This pattern of bacteria is associated with impaired intestinal barrier function (low fecal sIgA and EPX). Patients in this zone have higher rates of opportunistic infections (e.g. *Blastocystis* spp. & *Dientamoeba fragilis*) as well as fecal fat malabsorption. Commensal abundance is higher in this group suggesting potential bacterial overgrowth.

Zone 3: Patients in this zone may have more inflammation compared to those in zone 4. However, commensal abundance is usually higher making use of antimicrobial therapy relatively safer. Patients in this zone may have higher rates of pathogenic infections.


Zone 4: This commensal profile is associated with increased intestinal inflammation. IBD patients are more likely to have this pattern of bacteria. Commensal abundance is lower in this zone; therefore, antibiotic use for GI potential pathogens should be used with caution. In addition to standard treatment for intestinal inflammation, modulation of the commensal gut profile is encouraged.





Physician Notes/Recommendations

Commensal Abundance

- **Shift-to-the-Right:** Patient has more overall commensal bacteria
 - May be indicative of potential microbial overgrowth, such as in small intestinal bacterial overgrowth (SIBO)
 - May also be due to recent supplementation with probiotics
- **Shift-to-the-Left:** Patient has less overall commensal bacteria
 - May be indicative of potential microbiome deficiency, such as following antibiotic use
 - May indicate a diet low in fiber and prebiotic foods



Dysbiosis Patterns



Dysbiosis Patterns: Genova's data analysis has led to the development of unique dysbiosis patterns, related to key physiologic disruptions, such as immunosuppression and inflammation. These patterns may represent dysbiotic changes that could pose clinical significance. Please see Genova's published literature for more details: <https://rdcu.be/bRhzy>

Zone 1: The commensal profile in this zone does not align with profiles associated with intestinal inflammation or immunosuppression. If inflammatory biomarkers are present, other causes need to be excluded, such as infection, food allergy, or more serious pathology.

Zone 2: This pattern of bacteria is associated with impaired intestinal barrier function (low fecal sIgA and EPX). Patients in this zone have higher rates of opportunistic infections (e.g. *Blastocystis spp.* & *Dientamoeba fragilis*) as well as fecal fat malabsorption. Commensal abundance is higher in this group suggesting potential bacterial overgrowth.

Zone 3: Patients in this zone may have more inflammation compared to those in zone 4. However, commensal abundance is usually higher making use of antimicrobial therapy relatively safer. Patients in this zone may have higher rates of pathogenic infections.

Zone 4: This commensal profile is associated with increased intestinal inflammation. IBD patients are more likely to have this pattern of bacteria. Commensal abundance is lower in this zone; therefore, antibiotic use for GI potential pathogens should be used with caution. In addition to standard treatment for intestinal inflammation, modulation of the commensal gut profile is encouraged.


Inflammation-Associated Dysbiosis Score (IAD)

- Specific dysbiosis pattern associated with inflammation
- Correlated with inflammatory biomarkers
 - Calprotectin
 - Eosinophil Protein X
 - Secretory IgA
- Algorithm-derived from commensal bacteria analysis



Methane Dysbiosis Score

- Specific dysbiosis pattern associated with *methane* production
- Correlated with methane production on Genova SIBO tests
- Based both on commensal bacterial profile and stool biomarkers
- Developed an algorithm-derived score to predict higher methane production in the GI tract




Dysbiosis Patterns



Zone 1: The commensal profile in this zone does not align with profiles associated with intestinal inflammation or immunosuppression. If inflammatory biomarkers are present, other causes need to be excluded, such as infection, food allergy, or more serious pathology.

Commensal Microbiome Analysis

Commensal Balance




Balanced	Represents 95% of healthy individuals
Borderline	Represents 5% of healthy individuals
Imbalanced	Represents 60% of unhealthy individuals


*A progressive ranking scale based on a Genova proprietary algorithm that differentiates healthy and unhealthy commensal patterns.

**The total number of commensal bacteria (qPCR) that are out of balance for this individual on a scale of 0 to >12.

Relative Commensal Abundance










Relative Abundance: The relative abundance compares the quantity of each of 7 major bacterial phyla to a healthy cohort. This can indicate broader variances in the patient's gut microbiome profile. Certain interventions may promote or limit individual phyla when clinically appropriate. Please refer to Genova's Stool Testing Support Guide for more information on modulation of commensal bacteria through diet & nutrient interventions. ***Approximately 70% of the healthy cohort had below detectable levels of *Methanobrevibacter smithii*. Approximately 90% of the healthy cohort had below detectable levels of *Fusobacterium* spp.



2200 GI Effects™ Comprehensive Profile - Stool

Methodology: GC/MS, Automated Chemistry, EIA

2200 GI Effects™ Comprehensive Profile - Stool		QUINTILE DISTRIBUTION					Reference Range
	Result	1st	2nd	3rd	4th	5th	
Digestion and Absorption							
Pancreatic Elastase 1 †	158 L						>200 mcg/g
Products of Protein Breakdown (Total*) (Valerate, Isobutyrate, Isovalerate)	6.0						1.8-9.9 micromol/g
Fecal Fat (Total*)	19.5						3.2-38.6 mg/g
Triglycerides	1.1						0.3-2.8 mg/g
Long-Chain Fatty Acids	12.9						1.2-29.1 mg/g
Cholesterol	0.5						0.4-4.8 mg/g
Phospholipids	5.0						0.2-6.9 mg/g

Inflammation and Immunology							
Calprotectin †	145 H						<=50 mcg/g
Eosinophil Protein X (EPX) †	4.9 H		1.1		4.6		<=4.6 mcg/g
Fecal secretory IgA	206						<=885 mcg/g

Gastrointestinal Microbiome							
Metabolic							
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	81.3						>=23.3 micromol/g
n-Butyrate Concentration	18.1						>=3.6 micromol/g
n-Butyrate %	22.3						11.8-33.3 %
Acetate %	63.1						48.1-69.2 %
Propionate %	14.6						<=29.3 %
Beta-glucuronidase	2,297						368-6,266 U/g

2200 GI Effects™ Comprehensive Profile - Stool

Methodology: GC-FID, Automated Chemistry, EIA

Result

QUINTILE DISTRIBUTION
1st 2nd 3rd 4th 5th

Reference Range

Digestion and Absorption

Pancreatic Elastase 1 †

>500



>200 mcg/g

Digestion and Absorption

Pancreatic Elastase 1 †

>500



>200 mcg/g

Products of Protein Breakdown (Total*)
(Valerate, Isobutyrate, Isovalerate)

2.2



1.8-9.9 micromol/g

Fecal Fat (Total*)

6.6



3.2-38.6 mg/g

Triglycerides

0.7



0.3-2.8 mg/g

Long-Chain Fatty Acids

4.6



1.2-29.1 mg/g

Cholesterol

0.8



0.4-4.8 mg/g

Phospholipids

0.5



0.2-6.9 mg/g

n-Butyrate %

22.9



11.8-33.3 %

Acetate %

59.2



48.1-69.2 %

Propionate %

18.1



<=29.3 %

Beta-glucuronidase

1,547




368-6,266 U/g



Pancreatic Elastase 1

- A digestive enzyme secreted by the pancreas providing insight into pancreatic exocrine function
- Not affected by transit time, though profuse watery stool samples may falsely lower PE-1 due to dilution
- Not affected by digestive enzyme supplementation
- PE-1 correlates with the gold-standard secretin-cerulean test
- Low levels associated with chronic pancreatitis, gallstones, gastric bypass, Celiac disease, Diabetes, IBD, obesity

Products of Protein Breakdown






Products of Protein Breakdown

- Dietary protein not digested or absorbed effectively by the small intestine may be exposed to anaerobes in the colon which ferment them into

Products of Protein Breakdown

- Elevations may reflect poor digestion/absorption of protein
- Dietary intake can influence elevated findings based on higher protein volume
- Conversely, lower levels may indicate effective digestion/absorption of protein and/or lower dietary intake
- May reflect hypochlorhydria
- Check for SIBO
- Check for inflammation or infection

Fecal Fats





Fecal Fats

- Your patients should include fats as an essential component of their diets, elevated fecal fats imply poor digestion/absorption of those fats
- Triglycerides and cholesterol make up most of our dietary fat intake
- Triglycerides are broken down to form LCFAs
- Elevated fecal fats can be caused by:
 - Exocrine pancreatic insufficiency
 - Bile salt insufficiency
 - Use of PPIs and hypochlorhydria



Digestion and Absorption

Pancreatic Elastase 1 †	>500		>200 mcg/g
Products of Protein Breakdown (Total*) (Valerate, Isobutyrate, Isovalerate)	2.2		1.8-9.9 micromol/g
Fecal Fat (Total*)	6.6		3.2-38.6 mg/g
Triglycerides	0.7		0.3-2.8 mg/g
Long-Chain Fatty Acids	4.6		1.2-29.1 mg/g
Cholesterol	0.8		0.4-4.8 mg/g
Phospholipids	0.5		0.2-6.9 mg/g

Inflammation and Immunology

Calprotectin †	<16		<=50 mcg/g
Eosinophil Protein X (EPX)†	<DL		<=2.7 mcg/g
Fecal secretory IgA	683		<=2,040 mcg/mL

Gut Microbiome Metabolites

Metabolic


Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	29.3		>=23.3 micromol/g
n-Butyrate Concentration	6.7		>=3.6 micromol/g
n-Butyrate %	22.9		11.8-33.3 %
Acetate %	59.2		48.1-69.2 %
Propionate %	18.1		<=29.3 %
Beta-glucuronidase	1,547		368-6,266 U/g



Calprotectin

- Released from the intestinal mucosa into the stool in intestinal inflammation
- Fecal calprotectin is useful in differentiating IBD from IBS and monitoring IBD treatment
- It is **not** a cancer marker
- It is not a substitute for a scope, but can certainly direct the physician to the usefulness of scoping the patient
- Calprotectin 50-120 mcg/g can be caused by infection, hx of IBD, chronic NSAID or PPI use
- Calprotectin >120 refer to GI specialist to rule out IBD, malignancy

Eosinophil Protein X (EPX)








Eosinophil Protein X (EPX)

- Elevated with immune-mediated food hypersensitivity, atopic dermatitis and food allergies
- Inflammatory Bowel Disease (IBD)
- Certain parasitic infections
- Microscopic colitis (dx requires histological analysis)
- Can be elevated in children younger than 4 years old

Fecal Secretory IgA

Inflammation and Immunology					
Calprotectin †	<16		50	120	<=50 mcg/g
Eosinophil Protein X (EPX)†	<DL		0.5	2.7	<=2.7 mcg/g
Fecal secretory IgA	683		680	2040	<=2,040 mcg/mL



Fecal Secretory IgA

- Recognized as a first line of defense in protecting the intestinal epithelium from enteric pathogens
- Examples include Celiac disease, colon cancer, infections, IBS
- Treat root causes of immune upregulation/inflammation
- Assess for intestinal permeability
- Assess food antibody testing, consider use of an elimination diet
- Low sIgA may reflect a loss of GI immune response resiliency



Digestion and Absorption			
Pancreatic Elastase 1 †	>500		>200 mcg/g
Products of Protein Breakdown (Total*) (Valerate, Isobutyrate, Isovalerate)	2.2		1.8-9.9 micromol/g
Fecal Fat (Total*)	6.6		3.2-38.6 mg/g
Triglycerides	0.7		0.3-2.8 mg/g
Long-Chain Fatty Acids	4.6		1.2-29.1 mg/g
Cholesterol	0.8		0.4-4.8 mg/g
Phospholipids	0.5		0.2-6.9 mg/g
Inflammation and Immunology			
Calprotectin †	<16		<=50 mcg/g
Eosinophil Protein X (EPX)†	<DL		<=2.7 mcg/g
Fecal secretory IgA	683		<=2,040 mcg/mL
Gut Microbiome Metabolites			
Metabolic			
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	29.3		>=23.3 micromol/g
n-Butyrate Concentration	6.7		>=3.6 micromol/g
n-Butyrate %	22.9		11.8-33.3 %
Acetate %	59.2		48.1-69.2 %
Propionate %	18.1		<=29.3 %
Beta-glucuronidase	1,547		368-6,266 U/g



Short-Chain Fatty Acids







Gut Microbiome Metabolites				
Metabolic				
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	29.3		>=23.3 micromol/g	
n-Butyrate Concentration	6.7		>=3.6 micromol/g	
n-Butyrate %	22.9		11.8-33.3 %	
Acetate %	59.2		48.1-69.2 %	
Propionate %	18.1		<=29.3 %	
Beta-glucuronidase	1,547		368-6,266 U/g	



Short-Chain Fatty Acids

- Acetate, proprionate and butyrate are produced by bacterial fermentation of dietary fiber and resistant starch
- They act to maintain intestinal barrier function
- Provide fuel for colonocytes
- Support commensal bacteria
- Modulated anti-inflammatory and antimicrobial activities

Beta-glucuronidase

Gut Microbiome Metabolites				
Metabolic				
Short-Chain Fatty Acids (SCFA) (Total*) (Acetate, n-Butyrate, Propionate)	29.3			≥ 23.3 micromol/g
n-Butyrate Concentration	6.7			≥ 3.6 micromol/g
n-Butyrate %	22.9			11.8-33.3 %
Acetate %	59.2			48.1-69.2 %
Propionate %	18.1			≤ 29.3 %
Beta-glucuronidase	1,547			368-6,266 U/g



Beta-glucuronidase

- An enzyme produced by colonocytes and intestinal bacteria
- Can promote recirculation of various hormones and toxins that would have been eliminated
- Its action can therefore increase circulating estrogens
- Research suggests an association with increased risk of colorectal and breast cancer
- Elevation caused by dysbiosis and a Western diet high in red meat and protein
- Therapeutic considerations include probiotics, dietary fiber, Calcium-D-glucarate (found in oranges, apples, grapefruit and cruciferous vegetables)
- Low-calorie and vegetarian diets
- Konjac noodles are known to inhibit the action of the enzyme

Commensal Bacteria (PCR)

- Commensals are not inherently pathogenic
- Pattern analysis allows for better interpretation of findings
- PCR is quantitative
- Individual bacteria have unique clinical associations and importance to GI health
- PDF lists each individually with information on each

Gastrointestinal Microbiome (PCR)							
Commensal Bacteria (PCR)	Result CFU/g stool	QUINTILE DISTRIBUTION					Reference Range CFU/g stool
		1st	2nd	3rd	4th	5th	
Bacteroidetes Phylum							
<i>Bacteroides uniformis</i>	3.5E8						<=9.5E8
<i>Phocaeicola vulgatus</i>	2.8E8						<=8.3E8
<i>Barnesiella spp.</i>	3.6E7						3.0E6-2.9E8
<i>Odoribacter spp.</i>	<DL						<=9.5E7
<i>Prevotella spp.</i>	1.2E9						6.6E7-3.8E9
Firmicutes Phylum							
<i>Anaerotruncus colihominis/massiliensis</i>	1.6E7						<=2.0E7
<i>Butyrivibrio crossotus</i>	<DL						<=3.3E7
<i>Clostridium spp.</i>	<DL						<=1.5E7
<i>Coprococcus eutactus</i>	<DL						<=1.2E8
<i>Faecalibacterium prausnitzii</i>	2.4E8						1.1E6-1.1E9
<i>Lactobacillus spp.</i>	5.6E3						<=1.6E6
<i>Pseudoflavonifractor spp.</i>	1.4E6						1.3E4-2.9E7
<i>Roseburia spp.</i>	7.4E7						3.6E5-4.6E8
<i>Ruminococcus bromii</i>	4.6E8						<=1.5E9
<i>Veillonella spp.</i>	4.6E5						<=4.1E6
Actinobacteria Phylum							
<i>Bifidobacterium spp.</i>	5.0E7						4.6E5-2.6E8
<i>Bifidobacterium longum subsp. longum</i>	<DL						<=1.3E8
<i>Collinsella aerofaciens</i>	<DL						<=1.3E8
Proteobacteria Phylum							
<i>Desulfovibrio piger</i>	<DL						<=5.4E7
<i>Escherichia coli</i>	2.1E4						<=7.5E6
<i>Oxalobacter formigenes</i>	<DL						<=1.1E7
Euryarchaeota Phylum							
<i>Methanobrevibacter smithii</i>	<DL						<=2.0E7
Fusobacteria Phylum							
<i>Fusobacterium spp.</i>	<DL						<=1.8E5
Verrucomicrobia Phylum							
<i>Akkermansia muciniphila</i>	5.9E5						>=8.5E3

Commensal Bacteria Guide

Commensal Bacteria


The most current, literature-based information on human studies related to increased or decreased levels of the commensal bacteria is summarized in the following chart. Note that the findings in the literature may not be consistent with Genova's findings due to different methodologies, thus treatment efficacy may vary. Most therapeutic interventions do not work in isolation, meaning they do not exclusively only target that one organism. Genova has not conducted outcome studies on the impact of certain therapeutics on the microbiome markers. Clinician discretion is advised for appropriateness of therapeutics.

Under certain conditions, environmental factors may influence specific commensals to become pathobionts. Pathobionts are distinguished from true infectious agents; they are potential pathogens under certain conditions. It is unknown whether these organisms play a causative role in disease or are a consequence of a disease state. Literature is evolving regarding the definition of a pathobiont and the role of commensal bacteria.¹⁻³

Organism	Description	Increased Levels	Decreased Levels
<i>Bacteroides uniformis</i>	<p><i>Bacteroides uniformis</i> is a fiber-degrading bacteria. It colonizes the gut in early infancy and is promoted by breast feeding.⁴</p> <p>Thought to enhance the gut barrier through the production of butyrate and GABA.^{5,6} Also produces beta glucuronidase, degrades mucin, and produces folate.^{4,7,8}</p> <p>Studied in preclinical trials as a potential probiotic for use in inflammatory and metabolic disorders.⁹⁻¹¹ <i>B. uniformis</i> was found to be decreased in obese patients as compared to healthy or lean groups.^{12,13} It was higher in healthy controls as compared to patients with ulcerative colitis.¹⁴</p> <p>Enriched in healthy individuals versus colorectal cancer patients.¹⁵</p> <p>Associated with degradation of the isoflavone genistein, which then becomes less bioavailable to the human.¹⁶</p>	<p>In ten healthy males, the consumption of red wine polyphenols for 4 weeks significantly increased the amount of <i>Bacteroides uniformis</i> as well as other commensal bacteria species.¹⁷</p> <p>Higher levels of insoluble fiber are associated with higher levels of <i>B. uniformis</i>.¹⁸</p> <p>A more favorable metabolic risk profile in men on a healthy plant-based diet was seen with a certain microbial profile featuring increased <i>B. uniformis</i> and decreased <i>Prevotella copri</i>. The healthy diet was characterized by a higher intake of fiber, plant proteins, whole grains, fruits, vegetables, nuts, and legumes, and a lower intake of energy, animal proteins, refined grains, potatoes, sweets, animal fat, egg, dairy, and meats.¹⁹</p> <p>A small study (n=13) showed the presence of <i>B. uniformis</i> and other <i>Bacteroides</i> species in non-vegetarians, versus vegetarians.²⁰</p>	<p>Higher fiber intake from beans is associated with lower abundance of <i>B. uniformis</i>.²¹</p>
<i>Phocaeicola vulgatus</i>	<p>Generally considered a beneficial gut commensal, although is capable of attaching to and invading colonic epithelial cells and inducing pro-inflammatory cytokines.²²</p> <p>Produces beta-glucuronidase, succinate, lactate, acetate, formate, and propionate.^{23,24}</p>	<p>A high beef diet was associated with increases in <i>Bacteroides fragilis</i>, <i>B. vulgatus</i> and <i>Clostridium</i> spp. in 10 volunteers.²⁷</p>	<p>Decreased levels were found in 7-12-year olds who consumed oligofructose-enriched inulin (<i>BENEO's</i> prebiotic fiber <i>Synergy1</i>) for 16 weeks in a double-blind-controlled trial.²⁸</p>





Stool Culture

- Culture means it is living, which means sensitivities can be used for treatment protocol design
- Genova distinguishes pathogens, potential pathogens and non-pathogen findings
- Mycology is culture specific to viable yeast growth



Gastrointestinal Microbiome (Culture)

Human microflora is influenced by environmental factors and the competitive ecosystem of the organisms in the GI tract. Pathogenic significance should be based upon clinical symptoms.

Microbiology Legend			
NG	NP	PP	P
			
No Growth	Non-Pathogen	Potential Pathogen	Pathogen

Bacteriology (Culture)

Lactobacillus spp.

Escherichia coli

Bifidobacterium (Anaerobic Culture)

Additional Bacteria

Salmonella spp.

Shigella spp.

Klebsiella pneumoniae

Bacillus species

Enterococcus faecium

Mycology (Culture)

Candida krusei


Yeast, not *Candida albicans*

Additional Bacteria

Non-Pathogen: Organisms that fall under this category are those that constitute normal, commensal flora, or have not been recognized as etiological agents of disease.

Potential Pathogen: Organisms that fall under this category are considered potential or opportunistic pathogens when present in heavy growth.

Pathogen: The organisms that fall under this category have a well-recognized mechanism of pathogenicity in clinical literature and are considered significant regardless of the quantity that appears in the culture.





Bacteria Sensitivity

Prescriptive Agents

<i>Klebsiella pneumoniae</i>	R	I	S-DD	S	NI
Ampicillin	R				
Amox./Clavulanic Acid				S	
Cephalothin				S	
Ciprofloxacin				S	
Tetracycline				S	
Trimethoprim/Sulfa				S	

Natural Agents

<i>Klebsiella pneumoniae</i>	LOW INHIBITION	HIGH INHIBITION
Berberine		
Oregano		
Uva-Ursi		

Mycology Sensitivity

Candida Susceptibility Profile for Azoles*






Organism	Number of Isolates	% Sensitive	
		Fluconazole	Voriconazole
<i>Candida albicans</i>	25561	99.19%	99.51%
<i>Candida parapsilosis</i>	8777	98.64%	99.33%
<i>Candida krusei</i>	3420	0.23%	97.79%
<i>Candida tropicalis</i>	1076	93.22%	90.57%
<i>Candida glabrata</i>	2898	27.1%	90.9%

***Results of pharmaceutical sensitivities against certain yeast species are based on internal Genova data pertaining to the frequency of susceptibility of the specific yeast to the listed antifungal agent. The pharmaceutical results are not patient-specific. Conversely, the results of inhibition to nystatin and natural agents are patient-specific.**

Non-absorbed Antifungals

<i>Candida krusei</i>	LOW INHIBITION	HIGH INHIBITION
Nystatin		

Natural Agents

<i>Candida krusei</i>	LOW INHIBITION	HIGH INHIBITION
Berberine		
Caprylic Acid		
Garlic		
Undecylenic Acid		
Uva-Ursi		

Microscopic O&P Parasite

- Choice of one specimen or 3 over different days to cast a wider net
- Microscopic exam allows for a much wider capacity to identify parasites
- Includes WBC and Charcot-Leyden Crystals
- Findings of Few, Moderate, Many

Parasitology

Microscopic O&P Results

Microscopic O&P is capable of detecting all described gastrointestinal parasites. The organisms listed in the box represent those commonly found in microscopic stool analysis. Should an organism be detected that is not included in the list below, it will be reported in the Additional Results section. These results were obtained using wet preparation(s) and trichrome stained smear. For an extensive reference of all potentially detectable organisms, please visit www.gdx.net/product/gi-effects-comprehensive-stool-test

Genus/species	Result
Nematodes - roundworms	
<i>Ancylostoma/Necator</i> (Hookworm)	Not Detected
<i>Ascaris lumbricoides</i>	Not Detected
<i>Capillaria philippinensis</i>	Not Detected
<i>Enterobius vermicularis</i>	Not Detected
<i>Strongyloides stercoralis</i>	Not Detected
<i>Trichuris trichiura</i>	Not Detected
Cestodes - tapeworms	
<i>Diphyllobothrium latum</i>	Not Detected
<i>Dipylidium caninum</i>	Not Detected
<i>Hymenolepis diminuta</i>	Not Detected
<i>Hymenolepis nana</i>	Not Detected
<i>Taenia</i> spp.	Not Detected
Trematodes - flukes	
<i>Clonorchis/Opisthorchis</i> spp.	Not Detected
<i>Fasciola</i> spp./ <i>Fasciolopsis buski</i>	Not Detected
<i>Heterophyes/Metagonimus</i>	Not Detected
<i>Paragonimus</i> spp.	Not Detected
<i>Schistosoma</i> spp.	Not Detected
Protozoa	
<i>Balantidium coli</i>	Not Detected
<i>Blastocystis</i> spp.	Many Detected
<i>Chilomastix mesnili</i>	Not Detected
<i>Cryptosporidium</i> spp.	Not Detected
<i>Cyclospora cayentanensis</i>	Not Detected
<i>Dientamoeba fragilis</i>	Not Detected
<i>Entamoeba coli</i>	Not Detected
<i>Entamoeba histolytica/dispar</i>	Not Detected
<i>Entamoeba hartmanii</i>	Not Detected
<i>Entamoeba polecki</i>	Not Detected
<i>Endolimax nana</i>	Not Detected
<i>Giardia</i>	Not Detected
<i>Iodamoeba buetschlii</i>	Not Detected
<i>Cystoisospora</i> spp.	Not Detected
<i>Trichomonads</i> (e.g. <i>Pentatrichomonas</i>)	Not Detected
Additional Findings	
White Blood Cells	Not Detected
Charcot-Leyden Crystals	Not Detected
Other Infectious Findings	

PCR Parasitology

- PCR detection can only find those parasites the test is designed specifically to identify
- Combining PCR with Microscopic O&P provides a much wider possibility for detection
- Why don't we provide sensitivities for parasites?

Parasitology				
PCR Parasitology - Protozoa				Methodologies: DNA by PCR
Organism	Result	Units		Expected Result
<i>Blastocystis</i> spp.	<2.14e2	femtograms/microliter C&S stool	Detected	Not Detected
<i>Cryptosporidium parvum/hominis</i>	<1.76e2	genome copies/microliter C&S stool	Not Detected	Not Detected
<i>Cyclospora cayetanensis</i>	<2.65e2	genome copies/microliter C&S stool	Not Detected	Not Detected
<i>Dientamoeba fragilis</i>	<1.84e2	genome copies/microliter C&S stool	Not Detected	Not Detected
<i>Entamoeba histolytica</i>	<9.64e1	genome copies/microliter C&S stool	Not Detected	Not Detected
<i>Giardia</i>	<1.36e1	genome copies/microliter C&S stool	Not Detected	Not Detected

Additional Results		
Methodology: Fecal Immunochemical Testing (FIT)		
	Result	Expected Value
Fecal Occult Blood*	Negative	Negative
Color††	Brown	
Consistency††	Formed/Normal	

††Results provided from patient input.
 Tests were developed and their performance characteristics determined by Genova Diagnostics. Unless otherwise noted with *, the assays have not been cleared by the U.S. Food and Drug Administration.

Zonulin Family Peptide			
Methodology: EIA			
	Result	Reference Range	Zonulin Family Peptide
Zonulin Family Peptide, Stool	86.0	22.3-161.1 ng/mL	<p>This test is for research use only. Genova will not provide support on interpreting the test results. This test does not detect zonulin.¹ The Scheffler paper suggests that the IDK kit may detect a zonulin family peptide, such as properdin. Genova's unpublished data demonstrated that the current IDK kit results were associated with stool inflammation biomarkers and an inflammation-associated dysbiosis profile.</p> <p>The performance characteristics of Zonulin Family Peptide have been verified by Genova Diagnostics, Inc. The assay has not been cleared by the U.S. Food and Drug Administration.</p>

Additional Tests

Macroscopic/Direct Exam for Parasites

Methodology: Macroscopic Evaluation

No human parasite detected in sample.

Add-on Testing

Methodology: EIA

	Result	Expected Value
HpSA - <i>H. pylori</i>	Negative	Negative
<i>Campylobacter</i> spp.♦	Negative	Negative
<i>Clostridium difficile</i> ♦	Negative	Negative
Shiga toxin <i>E. coli</i> ♦	Negative	Negative
Fecal Lactoferrin♦	Negative	Negative




Key Points

- Evaluation of gut health is a key factor in all aspects of patient health
- Providing culture for bacteria and yeast detects *living* organisms, allowing for sensitivities to enhance protocol design
- Parasitology that includes a microscopic and PCR platform casts a much wider net for detection of parasites which shed unpredictably
- Genova's statistical analysis of Commensal Bacterial findings enhances pattern analysis and interpretation of findings



Support Materials



Commensal Bacteria

The most cur
or decreased
Note that the
due to differ
intervention:
that one orgi
certain thera
appropriater

Organ
Bacteroides-Prev

Pathogenic Bacteria & Yea

Genus/Organism

Aeromonas
Aeromonas hydrophila
Aeromonas caviae
Aeromonas veronii
Aeromonas jandaei
Aeromonas schuberti
Bacillus anthracis

Bacillus cereus

Parasitic Organ

NEMATODES – ROUNDWORMS

Organism

Ancylostoma -Necator
Ancylostoma duodenale
Necator americanus

Ascaris lumbricoides

Capillaria philippinensis

Enterobius vermicularis

Descriptio

Hookworms
Soil-transmitted
nematodes

(P)

Soil-transmitted
nematode
Most common hu
worm infection

(P)

(P)

Pinworm
The most common
infection in childn
5-10 in the US

(P)

Pathogen (P), Potential path

Pathogen (P), Potential pathogen (PP), Non-pathogen (N)

The New Advancing C

25
Mar
2020

Share This

