


GI Effects Comprehensive Stool Profile
Inflammation & Immunology: The Barrier to the Outside World

**Donna J Smith, ND, L.Ac**Medical Education Specialist | Genova Diagnostics





### **Functional Imbalance Scores**






## **GI Effects Biomarkers For Inflammation & Immunology**










## Calprotectin

- Released from the intestinal mucosa, indicates neutrophil mediated inflammation
- Useful in differentiating IBS from IBD, ex. Ulcerative Colitis and Crohn's, and monitoring IBD treatment
- Considerations when elevated:
  - o IBD, not in remission
  - Colorectal cancer and polyps
  - NSAID and PPI use
  - o Infection
  - o Diverticular disease
  - Bariatric surgery
  - Age, younger than 5 or greater than 60



# **Calprotectin**



Calprotectin (mcg/g stool)	Interpretation	Follow up
≤50	Normal range	None
50 to 120	Suggestive of low grade inflammation	Address cause of inflammation (Suspected or history of IBD, Infection, Chronic NSAID or PPI use, Diverticular disease) Consider anti-inflammatories such as fish oil Re-evaluate at 4-6 weeks
>120	Abnormal	Refer to GI specialist to rule out IBD, malignancy, or other cause of significant GI inflammation



## Next...

• Other markers to consider if Calprotectin is elevated.









## **Fecal Occult Blood**

- Blood not evident to the naked eye and present in microscopic quantities only
- Inflammation
- Infection
- Note: Bleeding hemorrhoids or menstruation, may interfere with testing.







## **White Blood Cells**

#### **Parasitology**

#### Microscopic O&P Results

Microscopic O&P is capable of detecting all described gastrointestinal parasites. The organisms listed in the box represent those commonly found in microscopic stool analysis. Should an organism be detected that is not included in the list below, it will be reported in the Additional Results section. These results were obtained using wet preparation(s) and trichrome stained smear. For an extensive reference of all potentially detectable organisms, please visit <a href="www.gdx.net/product/gi-effects-comprehensive-stool-test">www.gdx.net/product/gi-effects-comprehensive-stool-test</a>

Genus/species	Result	
Nematodes - roundworms		
Ancylostoma/Necator (Hookworm)	Not Detected	
Ascaris lumbricoides	Not Detected	
Capillaria philippinensis	Not Detected	
Enterobius vermicularis	Not Detected	
Strongyloides stercoralis	Not Detected	
Trichuris trichiura	Not Detected	
Cestodes - tapeworms		
Diphyllobothrium latum	Not Detected	
Dipylidium caninum	Not Detected	
Hymenolepis diminuta	Not Detected	
Hymenolepis nana	Not Detected	
Tagnia enn	Not Detected	

#### **Additional Findings**

White Blood Cells Charcot-Leyden Crystals Not Detected Not Detected

Balantidium coli	Not Detected
Blastocystis spp.	Many Detected
Chilomastix mesnili	Not Detected
Cryptosporidium spp.	Not Detected
Cyclospora cayetanensis	Not Detected
Dientamoeba fragilis	Not Detected
Entamoeba coli	Not Detected
Entamoeba histolytica/dispar	Not Detected
Entamoeba hartmanii	Not Detected
Entamoeba polecki	Not Detected
Endolimax nana	Not Detected
Giardia	Not Detected
Iodamoeba buetschlii	Not Detected
Cystoisospora spp.	Not Detected
Triphomonada (o.a. Pontotriphomonas)	Not Detected

#### **Additional Findings**

White Blood Cells Not Detected Charcot-Leyden Crystals Not Detected





# White Blood Cells (WBC)

Indicates an immune response that can be seen in

infectious conditions or inflammatory bowel disease (IBD).





## **Brief Review of Immunoglobulins**

There are 5 types of immunoglobulins, also known as **antibodies IgA, IgD, IgE, IgG, and IgM** 







#### ----

# **Immunoglobulins**

#### IgA (dimer)







### Immunoglobulin A

- First line defense against foreign pathogens
- Found in fluids made by mucus membranes
  - tears, saliva, respiratory and intestinal secretions, and colostrum.

#### Immunoglobulin D


- Blood in small amounts
- Function not clearly understood

### • Immunoglobulin E

- Associated with allergic reactions
- Stimulates histamine release



## **Immunoglobulins**



### Immunoglobulin G

- Protects against bacteria and viral infections
  - Secondary immune response
- Most abundant circulating immunoglobulin
- Found in blood and tissue fluids
- Immunoglobulin M
  - 1st antibody made to fight a new infection
  - Blood and lymph







- Immediate onset (minutes to hours)
- Symptoms of Allergic Reaction
   Hives, stuffy or itchy nose, sneezing, itchy, teary eyes, vomiting, stomach cramps or diarrhea, angioedema or swelling, shortness of breath or wheezing, anaphylaxis
- 90% of all food allergy reactions are cause by 8 foods

Milk, soy, eggs, wheat, peanuts, tree nuts, fish, and shellfish

Treatment suggestion
 Lifetime elimination
 Allergen-specific immunotherapy





- Delayed onset, hours to days
- Symptoms of food sensitivities
  - Headache, itching, bloating, fatigue, diarrhea, constipation, joint pain, indigestion, or worsening of your chronic health complaints
- Often associated with "Leaky Gut"
- Treatment suggestion
  - Elimination Diet
  - Intestinal Permeability
  - Digestion/Absorption
  - Dysbiosis
  - Nutrition






# **Eosinophil Protein X (EPX)**

Inflammation and Immunology				
Calprotectin †	<16	50	120	<=50 mcg/g
Eosinophil Protein X (EPX)†	<dl< td=""><td>0.5</td><td>2.7</td><td>&lt;=2.7 mcg/g</td></dl<>	0.5	2.7	<=2.7 mcg/g
Fecal secretory IgA	683	680 ◆	2040	<=2,040 mcg/n



# **Eosinophil Protein X (EPX)**




- Eosinophil mediated inflammation
- Elevated with IgE mediated food allergies, atopic dermatitis, and immune-mediated food hypersensitivity
- Inflammatory Bowel Disease (IBD)
- Certain parasitic infections
- Microscopic colitis/gastritis (dx requires histological analysis)
- Eosinophilic gastrointestinal disorders
- Can be elevated in children younger than 4 years old





- IgE-mediated allergy
   IgE Food Antibody panel
- IBD (review Calprotectin level)
- Evaluate for parasitic infection
- Anti-Inflammatories
- Mast Cell Stabilizers





# **Fecal Secretory IgA**

Calprotectin †	<16	<b>5</b> 0	120	<=50 mcg/g
Eosinophil Protein X (EPX)†	<dl< td=""><td>0.5</td><td>2.7</td><td>&lt;=2.7 mcg/g</td></dl<>	0.5	2.7	<=2.7 mcg/g
Fecal secretory IgA	683	680	2040	<=2,040 mcg/mL



# Fecal Secretory IgA (fsIgA)




- Non-specific indicator of an immune response at the gut barrier
- Recognized as a first line of defense in protecting the intestinal epithelium from enteric pathogens
- Considerations when elevated:
  - Infection/Potential Pathogen
  - IgG Food Sensitivity, Celiac disease
  - Defective Epithelial Barrier, Gut Permeability
  - Colon cancer, IBS and dysbiosis
- Low slgA may reflect a loss of GI immune response resiliency



## Therapeutic Considerations for Elevated fslgA


- Address infection
- Assess for compromised intestinal barrier function
- Elimination Diet





## Therapeutic Considerations for Elevated fslgA

- IgG Food antibody testing
- Celiac testing






## Next...


Optional Add On Markers



## **Zonulin Family Peptide**

- Non-specific indicator associated with intestinal inflammation
- Clinical significance of an elevated Zonulin Family
   Peptide is unknown
- May relate to increased intestinal permeability.
   Consider follow up testing with the Intestinal Permeability Assessment.











### **Fecal Lactoferrin**

- Biomarker of neutrophil mediated inflammation.
- Expressed a Positive or Negative finding.
- Calprotectin is a more sensitive biomarker.






### **Functional Imbalance Scores**






## **Case Study 1**

- 39 yo M
- Hx of IBD,
   Ulcerative Colitis
- Carnivore diet




Detected
Not Detected





## Case Study 2

- 84 yo F
- Diarrhea
- Abdominal pain
- Weight loss
- Hives







- 64 yo M
- Urgent diarrhea
- Gas, cramping
- Recently returned from a 2 week backpacking trip

	In	nflammation and Immunology	
Calprotect	commonly found in microscopic stool analysis. Sh	ribed gastrointestinal parasites. The organisms listed in the box represent those nould an organism be detected that is not included in the list below, it will be reported	cg/g
Eosinophil	in the Additional Results section. These results w reference of all potentially detectable organisms,	ere obtained using wet preparation(s) and trichrome stained smear. For an extensive please visit www.gdx.net/product/gi-effects-comprehensive-stool-test	cg/g
Fecal secr	Genus/species	Result	) mcg/m
	Nematodes - roundworms		9
	Ancylostoma/Necator (Hookworm)	Not Detected	
	Ascaris lumbricoides	Not Detected	
	Capillaria philippinensis	Not Detected	
	Enterobius vermicularis	Not Detected	
	Strongyloides stercoralis	Not Detected	
	Trichuris trichiura	Not Detected	
	Cestodes - tapeworms		
	Diphyllobothrium latum	Not Detected	
	Dipylidium caninum	Not Detected	
	Hymenolepis diminuta	Not Detected	


PCR Parasitology - Proto	Methodologies: DNA by PCR			
Organism	Result	Units		Expected Result
Blastocystis spp.	<2.14e2	femtograms/microliter C&S stool	Not Detected	Not Detected
Cryptosporidium parvum/hominis	<1.76e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Cyclospora cayetanensis	<2.65e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Dientamoeba fragilis	<1.84e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Entamoeba histolytica	<9.64e1	genome copies/microliter C&S stool	Not Detected	Not Detected
Giardia	9.48e2	genome copies/microliter C&S stool	Detected	Not Detected


Not Detected
Not Detected
Not Detected
Rare Cyst(s) Detected
Not Detected
Not Detected
Not Detected
Not Detected
Not Detected



## Case Study 4

- 27 yo F
- Alternating constipation, diarrhea
- Gas, bloating
- Preparing for pregnancy
- Takes PERT







## Case Study 4 cont.

- 27 yo F
- Alternating constipation, diarrhea
- Gas, bloating
- Preparing for pregnancy
- Takes PERT

Parasitology					
Microscopic O&P Results  Microscopic O&P is capable of detecting all described gastrointestinal parasites. The organisms listed in the box represent those commonly found in microscopic stool analysis. Should an organism be detected that is not included in the list below, it will be reported in the Additional Results section. These results were obtained using wet preparation(s) and trichrome stained smear. For an extensive reference of all potentially detectable organisms, please visit  www.gdx.net/product/gi-effects-comprehensive-stool-test					
Genus/species	Result				
Nematodes - roundworms					
Ancylostoma/Necator (Hookworm)	Not Detected				
Ascaris lumbricoides	Not Detected				


PCR Parasitology - Proto	Methodologies: DNA by PCR			
Organism	Result	Units		Expected Result
Blastocystis spp.	1.99e4	femtograms/microliter C&S stool	Detected	Not Detected
Cryptosporidium parvum/hominis	<1.76e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Cyclospora cayetanensis	<2.65e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Dientamoeba fragilis	<1.84e2	genome copies/microliter C&S stool	Not Detected	Not Detected
Entamoeba histolytica	<9.64e1	genome copies/microliter C&S stool	Not Detected	Not Detected
Giardia	<1.36e1	genome copies/microliter C&S stool	Not Detected	Not Detected

Protozoa	
Balantidium coli	Not Detected
Blastocystis spp.	Many Detected
Chilomastix mesnili	Not Detected
Cryptosporidium spp.	Not Detected
Cyclospora cayetanensis	Not Detected
Dientamoeba fragilis	Not Detected
Entamoeba coli	Not Detected
Entamoeba histolytica/dispar	Not Detected
Entamoeba hartmanii	Not Detected
Entamoeba polecki	Not Detected
Endolimax nana	Not Detected
Giardia	Not Detected
Iodamoeba buetschlii	Not Detected
Cystoisospora spp.	Not Detected
Trichomonads (e.g. Pentatrichomonas)	Not Detected



## Case Study 4 cont.

- 27 yo F
- Alternating constipation, diarrhea
- Gas, bloating
- Preparing for pregnancy
- Takes PERT







- 27 yo F
- Alternating constipation, diarrhea
- Gas, bloating
- Preparing for pregnancy
- Takes PERT

Inflammation and Immunology				
Calprotectin †	<16	50	120	<=50 mcg/g
Eosinophil Protein X (EPX)†	<dl< td=""><td>0.5</td><td>2.7</td><td>&lt;=2.7 mcg/g</td></dl<>	0.5	2.7	<=2.7 mcg/g
Fecal secretory IgA	<150	680	2040	<=2,040 mcg/mL



### **Additional Educational Resources**

#### www.GDX.net


- GI Effects Support Guide and organism charts
- Learning Library video modules
- Test Prep information

#### **The Lab Report Podcast**

Available on Apple Podcasts and GDX.net

#### **Medical Education Consultations**

- Schedule online through myGDX
- Call Client Services 800-522-4762
- 1:1 and group consultations









GI Effects Comprehensive Stool Profile
Inflammation & Immunology: The Barrier to the Outside World

**Donna J Smith, ND, L.Ac**Medical Education Specialist | Genova Diagnostics

